

Active®

The Outdoor Rollator

The Active is the ultimate outdoor rollator

Sturdy and comfortable, with large front wheels, the Active takes you anywhere you want to go. With a comfortable seat to sit down on when you need a break, a convenient, secure shopping basket and, of course, the safety features you need to be able to enjoy an active life both indoors and out.



Safety & comfort - made easy

A rollator should be a reliable companion, and we have designed a product that is safe and comfortable to use in any daily situation, for just about anyone.



The Active was designed in Norway, by an Access team dedicated to safety, comfort and the possibility of maintaining an active and independent lifestyle.



Large wheels – Perfect for outdoors and most terrains, keeping the Active stable and safe, even on rough ground.

Snag-free cables – The brake wires are hidden inside the frame to keep from snagging, to keep you safe and to avoid dragging items from surfaces as you pass by.



Folding – One handed and easy. The Active stands folded so you can park it anywhere you want.

Fits all – Active has one of the widest ranges of height adjustments, to make sure you find the correct height for you.



Kerb climber – We recommend using our convenient kerb climber which saves you lifting the rollator up onto kerbs and over thresholds.

Room to move – Due to the ergonomic design of the Active, you have plenty of room for long, upright, comfortable strides. Many other rollators encourage you to lean forward, which is bad for your back.





Easy access - To get close to counters etc, just lift up the seat, remove the basket and you are inside the rollator with support and comfort, as with every other situation using the Active.

Seat height – Because we are all different, the seat is adjustable, to make it comfortable and ergonomically correct. And of course it is as easy to adjust as everything else on the Active.



"Safety in numbers" - Our height levels are numbered and click into place so that you can be sure to have the correct height, and that the position is locked and safe.

Easy adjustments – Everything is easy to adjust! No tools needed and no loose parts that can fall away and get lost. Every adjustment is marked in red and is intuitive and easy to do.



Safe storage - We have placed the shopping basket under the seat – to keep your belongings safe. *The Active comes fully assembled, so all you need to do is take it out of the box and you are ready to go!*

“I ordered the Active Rollator prior to having bilateral knee replacement surgery. It has been a wonderful tool in my recovery”

Williams, US



AWARD FOR DESIGN EXCELLENCE
NORWEGIAN DESIGN COUNCIL

Did you know that...

The Active was awarded the Norwegian award for design excellence.

Here is the verdict:

“The Active rollator has solved many challenges in a better way than that of its competitors”... “Fine attention to detail”... “a well thought through product”... “This product has been carefully designed with a sporty appearance adapted to a demanding and active user group”

We hope you will try the Active yourself and agree!

There is an International Standard for rollators (ISO 11199-2) requiring that “rollators manufactured for outdoor use shall have front wheel diameter no less than 180 mm”.

Going beyond the standards, the Active features 250 mm front wheels – for your safety.

The Active has a large range of accessories. Here are just some:

Active Cushion Set

Personalise your Active with these bright and colourful cushion sets. Available in seven different colours!



Active Deluxebag

Conveniently located in front of the seat, with easy attachment and a practical handle, this Deluxebag is perfect for increased storage space.



Active Backrest

For that extra support, this practical backrest moves with the seat to make sure you are always comfortable.



www.active-walker.com

Visit our international website for more information on the Active, including demonstration videos that will show just how easy it is to use and where to buy!

ACCESS
VITAL

T: +47 405 66 33 6

E: rollator@accessvital.no

W: active-walker.com

Technical specifications:

Type	Access Active
Dimensions	60 x 72 x 78-100cm
Folded dimensions	60 x 46 x 85cm
Weight	9
User weight	125kg
Front wheel size	25cm
Back wheel size	20cm
Seat height	50-63cm

Don't sit inside and watch the world go by – get Active and stay active!



Active[®]
by **ACCESS**